

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

- 1) 6 mmol/L of urea is equivalent to:-
 - a) **360mg/L**
 - b) 60mg/L
 - c) 36mg/L
 - d) 360g%

- 2) The product of protein deamination is:-
 - a) Creatinine
 - b) Creatine
 - c) **Used in clearance tests**
 - d) Creatine anhydride

- 3) Parameters used to assess liver function include:-
 - a) **Albumin level**
 - b) Plasma urea
 - c) Lipase activity
 - d) Trypsin activity

- 4) Ketoacidosis is associated with:-
 - a) Diabetes insipidus
 - b) Raised plasma pH
 - c) Post-prandial hyperglycaemia
 - d) **Acute starvation**

- 5) False positive results with albusix may be due to:-
 - a) **Alkaline urine**
 - b) Albuminuria
 - c) Acidic urine
 - d) High altitude

- 6) Detection of CSF total proteins is achieved by:-
 - a) Albustix
 - b) **Turbidimetric methods**
 - c) Using phenol solutions
 - d) Pandy test

- 7) In gastric juice analysis, free acidity refers to:-
 - a) Combined acid
 - b) Strong acid

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

- c) Total acid
 - d) HCl**
- 8) The highest level a substance is present in plasma before it is excreted is:-
- a) Indicative of diabetes mellitus
 - b) Indicative of diabetes insipidus
 - c) Renal threshold value**
 - d) Known as glucosuria
- 9) Biuret test involves the reaction of copper ions with:-
- a) 2-OH
 - b) 2C=O
 - c) 2-CONH**
 - d) 2-NH₂
- 10) Which of these can be used to make equivalent standards of urea:-
- a) Sodium sulphate
 - b) Potassium sulphate
 - c) Ammonium chloride**
 - d) Barium sulphate
- 11) Decarboxylated diacetic acid is known as:-
- a) Aceto acetic acid
 - b) Principal ketone body
 - c) Beta-hydroxy butyric acid
 - d) Acetone**
- 12) The separation phases in paper chromatography are:-
- a) Normal phase
 - b) Fast phase
 - c) Stationary phase**
 - d) Second phase
- 13) Occult test detects:-
- a) The presence of WBC's in urine
 - b) Malignancies of gastro-intestinal tract**
 - c) Peroxidase activity
 - d) Presence of dyes in stool

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

- 14) Which of the following plasma substance(s) will decrease in liver failure:-
a) Glucose
b) Creatinine
c) Uric acid
d) Transaminase
- 15) Which of the following is/are pre-analytical error(s):-
a) Reagent quality
b) Specimen collection technique
c) OD reading
d) Reporting of test result
- 16) Simple lipids are derived from:-
a) Esters
b) Alcohol
c) Acetone
d) Carboxylic acids
- 17) The following vitamin is required in blood clotting:-
a) B₆
b) C
c) B₁₂
d) K
- 18) Water prepared by distillation is:-
a) Sterile but has soluble salts
b) Sterile and free from soluble salts
c) Not sterile but free from soluble salts
d) Sterile only
- 19) To prepare 10ml of a 1:200 serum dilution, you need:-
a) 2ml of serum
b) 0.05ml of serum
c) 1 ml of serum
d) 0.5 ml of water
- 20) Glycosuria with hyperglycaemia denotes:-
a) Diabetes mellitus
b) Diabetes insipidus

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

- c) Lowered renal threshold
 - d) Renal glycosuria
- 21) Tetrabromophenol blue is an indicator in:-
- a) Clinistix
 - b) Urostrat
 - c) Glucostix
 - d) Albustix**
- 22) At pH below their isoelectric points, proteins:-
- a) Ionize as bases
 - b) Acquire zwitterionic state
 - c) Acquire cationic state**
 - d) Are anions
- 23) The direct reacting bilirubin is also:-
- a) Unconjugated
 - b) Insoluble in water
 - c) Water soluble**
 - d) Very toxic
- 24) In indirect agglutination test, agglutination is interpreted as:-
- a) Positive
 - b) False positive
 - c) False negative
 - d) Negative**
- 25) Surface tension of urine can be lowered by:-
- a) Bile salts**
 - b) Glucose
 - c) Amino acids
 - d) Bilirubin
- 26) What is the waste end product of purine metabolism:-
- a) Urea
 - b) Uric acid**
 - c) Creatine
 - d) Creatinine
- 27) Characteristics of enzymes include:-

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

- a) Activated by fluoride
 - b) Reacting at 37°C
 - c) Specificity**
 - d) Resists UV denaturation
- 28) Blood glucose levels are directly regulated by:-
- a) ADH
 - b) Cortisol
 - c) Glucagon**
 - d) Epinephrine
- 29) Oxidative deamination of proteins yields:-
- a) Urea**
 - b) Amino acids
 - c) Peptide residues
 - d) Peptide bonds
- 30) Amylase activity assay are performed to diagnose:-
- a) Gastric disorders
 - b) Diabetes mellitus
 - c) Chronic pancreatitis
 - d) Acute pancreatitis**
- 31) A control in clinical chemistry is:-
- a) Same as a standard solution
 - b) A solution of known reference range**
 - c) A solution of known dilution
 - d) A solution used to calibrate a machine
- 32) Parenteral stimulants of gastric juice secretion include:-
- a) Alcohol
 - b) Gruel
 - c) Insulin**
 - d) Tea without milk
- 33) The ability of plasma proteins to buffer blood is due to their ability to:-
- a) Neutralize H^+ and OH^-**
 - b) Chelate Ca^{2+}
 - c) Form peptide bonds
 - d) Form hydrogen bonds

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

- 34) Markers of cardiac diseases is/are:-
a) Lipase assay
b) Alkaline phosphatase activity assay
c) Raised tryptic activity values
d) Raised creatine kinase activity
- 35) Major fluid compartments of the body include:-
a) Intracellular
b) Interstitial space
c) Extravascular
d) Intravascular
- 36) Vitamin intermediates with hormone-like activity include:-
a) B₁
b) E
c) K
d) D
- 37) Man and other primates excrete nitrogen in form of:-
a) Uric acid
b) Ammonia
c) Urea
d) Creatinine
- 38) Diabetes insipidus is characterized by:-
a) Urine of high specific gravity
b) Lack of antidiuretic hormone (ADH)
c) Excess insulin secretion
d) Deficiency of aldosterone
- 39) Select from the list below cation(s) that are tested for in serum:-
a) Sodium
b) Calcium
c) Caesium
d) Magnesium
- 40) Enzymes of clinical importance analysed in stool include:-
a) Amylase
b) Trypsin

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

- c) Lipase
- d) Diastase

41) Chloride shift in the body ensures:-

- a) Fluid balance
- b) Sodium balance
- c) Cationic/anionic balance**
- d) Potassium balance

1) Indicate the organ assessed by assaying each of the analytes listed below:-

- i) Glucose
- ii) Acid phosphatase
- iii) Urea
- iv) LDH₄
- v) PSA

(10

marks)

2) Match the substances listed in column A with their respective organ/tissue source in column

	A	B
1	CK _{BB}	Corpus luteum
2	Amylase	Adrenal medulla
3	Progesterone	Brain
4	Epinephrine	Pancreas
5	Insulin	

(10

marks)

3) Discuss the formation of ketone bodies.

(10 marks)

4) List two analytes under the following function test/profile and the unit of measure for each.

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

- i) Liver
 - ii) Lipid
 - iii) Heart
 - iv) Pancreas
 - v) Kidney
- marks)**

(10

5) State the principle behind the following serum estimation methods:-

- a. Uricase
- b. Biuret test

**(10
marks)**

6) State the principle behind the following serum estimation methods:-

- a. Sackett's
- b. Jaffe's reaction

**(10
marks)**

Answers

1. Indicate the organ assessed by blood assay of each of the analytes listed below:-

i) Glucose	Pancreas	
ii) Troponin I	Heart	
iii) Urea	Kidney	
iv) LDH ₄	Liver	
v) PSA	Prostate	(10 marks)

2. Match the substances listed in column A with their respective organ/tissue source in column B.

	A	B
1	CK _{BB}	Brain
2	Amylase	Pancreas
3	Progesterone	Corpus luteum
4	Epinephrine	Adrenal medulla
5	Insulin	Pancreas

(10 marks)

marks)

3. List two analytes under the following function test/profile and the unit of measure for each.

- Liver **ALT, AST, ALP, u/l**
- Lipid **HDL, LDL, VLDL, CHOLESTEROL, mg/dl or Mmol/L**
- Heart **Troponin I, Troponin T, Myoglobin, Creatinine Kinase, ng/ml**
- Pancreas **GLUCOSE, INSULIN, LIPASE, AMYLASE, mmol/l U/L**
- Kidney **urea, creatinine, Na, K, Cl, Mmol/l**

(10 marks)

4. Discuss the formation of ketone bodies. **(10 marks)**

Ketone bodies are produced mainly in the mitochondria of liver cells, and synthesis can occur in response to an unavailability of blood glucose, such as during fasting

Ketogenesis takes place in the setting of low glucose levels in the blood, after exhaustion of other cellular carbohydrate stores, such as glycogen. It can also take place when there is insufficient insulin e.g. in diabetes. The production of ketone bodies is then initiated to make available energy that is stored as fatty acids.

Fatty acids are enzymatically broken down in β -oxidation to form acetyl-CoA. Under normal conditions, acetyl-CoA is further oxidized by the citric acid cycle (TCA/Krebs cycle) and then by the mitochondrial electron transport chain to release energy.

However, if the amounts of acetyl-CoA generated in fatty-acid β -oxidation challenge the processing capacity of the TCA cycle; i.e. if activity in TCA cycle is low due to low amounts of intermediates such as oxaloacetate, acetyl-CoA is then used instead in biosynthesis of ketone bodies via acetoacetyl-CoA and β -hydroxy- β -methylglutaryl-CoA (HMG-CoA).

Steps in ketogenesis:

1. Catalyzed by enzyme thiolase, 2 molecules of acetyl-CoA condense to form acetoacetyl-CoA.
2. Catalyzed by enzyme HMG-CoA synthase, acetoacetyl-CoA combines with another molecule of acetyl-CoA to produce HMG-CoA (β -hydroxy β -methyl glutaryl-CoA).
3. HMG-CoA lyase cleaves HMG-CoA to produce acetoacetate and acetyl-CoA.
4. Acetoacetate can undergo spontaneous decarboxylation to form acetone.
5. Acetoacetate can be reduced by a dehydrogenase to β -hydroxybutyrate.

5.

i. Uricase

Used for serum uric acid estimation. Uricase enzyme digests uric acid to form allantoin, hydrogen peroxide and carbon dioxide. The decrease of absorption at 293nm is monitored and is related to the uric acid present.

ii. Biuret reaction

Used for serum protein estimation. Peptide bonds in protein molecular react with copper (ii) ions in alkaline medium to form a violet/blue coloured complex. The intensity of the colour is proportional to the amount of protein present.

6.

i. Jaffe's reaction

Used for serum creatinine estimation. Creatinine reacts with picric acid in alkaline medium to form a red colour. The intensity of the red colour is proportional to the amount of creatinine present

ii. Sackett's method

Used for serum cholesterol estimation. Cholesterol reacts with acetic anhydride -sulphuric acid mixture to form a green colour. The intensity of the green colour is proportional to the amount of cholesterol present

Practical

You are provided with a urine specimen labelled “K”. Use the reagents and apparatus provided to detect abnormal substances in the specimen.

- (i) Record your results
- (ii) Outline the principles of the positive tests performed.
- (iii) Comment on the results
- (iv) Suggest further tests

(30 Marks)

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

PRACTICAL PREPARATION (02) CHEMISTRY

Provide the following each candidate

- Urine specimen "K" containing
 - Glucose
 - Protein
 - salicylates
- Combur strips

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

Results

Show all the parameters

Parameter	Result	Marks
Glucose	Pos	1
Protein	Pos	1
Ketones	Pos	1
Others		1
		1
		1
		1
		1
		1
		1

- Principles (2)
 - G.O. D
 - albusix (2)
 - Gerhardt's (2)
- Comment (1)
 - Specimen contains Glucose (1)
 - Protein (1)
 - Salicylates/ acetoacetic acid (1)
 - The specimen could be from a patient with
 - (i) Diabetes mellitus (1)
 - (ii) Kidney problem (1)

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

- (iii) Aspirin poisoning (1)
- Tests -
 - o OGTT/ FBS for Diabetes (2)
 - o Boiling to differentiate salicylates and acetoacetate (2)
- Neatness (2)
- - Work evidence - (2)

TOTAL 30

points

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

Q.2. You are provided with a 200 mg/dl protein standard. Using the reagents and apparatus provided, prepare a C.S.F. Protein calibration curve to represent 0, 40, 80, 120, 160 and 200 mg/dl.

- a) Use the curve to read the concentration of C.S.F. sample 'S'
- b) Comment on the result

PROCEDURE

- a) Mix 3 ml of 3% TCA with 1 ml of the protein standards
- b) Mix 3 ml of 3% TCA with 1 ml of C.S.F.
- c) Stand the tubes at room temperature for 10 minutes
- d) Set up a blank containing 3 ml of 3% TCA and 1 ml distilled water
- e) Read the absorbances (O.D) of the test and standards against blank using a blue filter or at 450 nm wavelength.

Prep guide

- Distilled water 5 ml each
- Normal saline - 10 ml each
- Sample "S" (each) 5 ml of 50 mg % protein solution
- Rack 1 each
- Test tubes 7 each
- 1 ml pipette 2 each
- 5 ml pipette 1 each
- 200 mg% protein solution - 10ml each
- Graph paper
- Colorimeter with filter - 450 nm

Marking Guide

1. Use of the formula $RV \times RC/OC$ **2 marks**
2. Calculation of dilutions 0, 40, 80, 120, 160, 200 mg/dl **1 mark each = 6 marks**
3. Table of dilutions

Tube	0	4	8	12	16	20	S
		0	0	0	0	0	
Standards (ml)	1	1	1	1	1	1	-
Test (ml)	-	-	-	-	-	-	1
3% TCA (ml)	3	3	3	3	3	3	3
Total (ml)	4	4	4	4	4	4	4
OD							

2 marks

4. Plot a curve
 - a. Label X and Y axes **1 mark**
 - b. Title **1 marks**
 - c. Date **1 mark**
 - d. Machine used **1 mark**
 - e. Wavelength and light path **1 mark**
 - f. Prepared by **1 mark**
 - g. Table of OD and Concentration **2marks**
 - h. 4 points in line **2 marks**
 - i. Reading of "S" **50 mg/dl +/- 2 mg/dl 2 marks**
 - j. Neatness **2 marks**
 - k. Evidence **2 marks**
5. Comment
 - i. Reference range 15-45 mg/dl **1 marks**

DIPLOMA CLINICAL CHEM MAY 2016 MODERATED

ii. Above normal **1**
mark

iii. Conditions meningitis, intra cranial bleeding
2 marks

30 marks